

Anti-phospho-PKC zeta (Thr410)/PKC lambda (Thr412) Rabbit antibody



Catalog Number: 167962

Product name

Anti-phospho-PKC zeta (Thr410)/PKC lambda (Thr412) Rabbit antibody

Ig Type

IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

phospho-PKC zeta (Thr410)/PKC lambda (Thr412)
Rabbit monoclonal antibody

Dilutions

WB: 1/1000

IHC: 1/20

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr410 of human PKC zeta

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of PKC zeta (phospho-Thr410)/PKC lambda (phospho-Thr412) in HeLa, A549, HL-60 using PKC zeta (phospho-Thr410)/PKC lambda (phospho-Thr412) Rabbit mAb(1:1000 diluted)