

Anti-RING2/RING1B/RNF2 Rabbit antibody



Catalog Number: 168022

Product name

IgG

Anti-RING2/RING1B/RNF2 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/500-1/2000

RING2/RING1B/RNF2 Rabbit monoclonal antibody

IHC: 1/50-1/200

Preparation

ICC/IF: 1/50-1/200

Antigen: A synthesized peptide derived from human RING2 / RING1B / RNF2

FC: 1/50

Validations

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of RING2/RING1B/RNF2 in HeLa,A549,HL-60,U2OS,U87-MG cell lysates using RING2/RING1B/RNF2 Rabbit mAb(1:1000 diluted).Observed band size:38kDa.