

Anti-FE65 Rabbit antibody



Catalog Number: 168029

Product name

Anti-FE65 Rabbit antibody

Specificity

Human

Antibody description

FE65 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FE65

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications


WB, IHC-P, ICC/IF

Dilutions

WB: 1/1000

IHC: 1/201/20

Validations



Western blot detection of FE65 in A549,U251 cell lysates using FE65 Rabbit mAb(1:1000 diluted).Predicted band size:77kDa.Observed band size:100kDa.