

Anti-Cytokeratin 17 Rabbit antibody



Catalog Number: 168036

Product name

Anti-Cytokeratin 17 Rabbit antibody

Specificity

Human

Antibody description

Cytokeratin 17 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cytokeratin 17

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/2000

IHC: 1/20

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of Cytokeratin 17 in HeLa cell lysates using Cytokeratin 17 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:48kDa.