

# Anti-CaMKII alpha (5D9) Mouse antibody



Catalog Number: 168084

## Product name

Anti-CaMKII alpha (5D9) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

CaMKII alpha (5D9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CaMKII- $\alpha$  protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of CaMKII- $\alpha$  in Mouse Brain lysates using CaMKII- $\alpha$  mouse mAb(dilution 1:1000).Predicted band size:50kDa.Observed band size:50kDa.