

Anti-HARS (4H5) Mouse antibody



Catalog Number: 168088

Product name

Anti-HARS (4H5) Mouse antibody

Specificity

Human, Monkey

Antibody description

HARS (4H5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human HARS protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/2000

Validations



Western blot detection of HARS in T47D,3T3,COS7,C6 and HeLa cell lysates using HARS mouse mAb(dilution 1:2000).Predicted band size:56kDa.Observed band size:56kDa.