

# Anti-Aldose Reductase/AKR1B1 (4G9) Mouse antibody



Catalog Number: 168089

## Product name

Monoclonal

Anti-Aldose Reductase/AKR1B1 (4G9) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

Aldose Reductase/AKR1B1 (4G9) Mouse monoclonal antibody

## Dilutions

WB: 1/1000


## Preparation

Antigen: Purified recombinant human AKR1B1 protein fragments expressed in E.coli.

## Validations

## Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of AKR1B1 in C6,3T3,Jurkat and Hela cell lysates using AKR1B1 mouse mAb (1:1000 diluted). Predicted band size:36KDa. Observed band size:36KDa.

## Clonality