

# Anti-FEN1 (7H8) Mouse antibody



Catalog Number: 168116

## Product name

Anti-FEN1 (7H8) Mouse antibody

## Dilutions

WB: 1/1000

## Specificity

;ICC: 1/400

Human, Mouse, Rat, Monkey, Hamster

## Validations

## Antibody description

FEN1 (7H8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human FEN-1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


IgG1

## Applications

WB, ICC/IF



Western blot detection of FEN-1 in HeLa, Jurkat, 3T3, COS7, PC-12, C6, Raw246.7, Raji and CHO-K1 cell lysates using FEN-1 mouse mAb (1:1000 diluted). Predicted band size: 45KDa. Observed band size: 45KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using FEN-1 mouse mAb (dilution 1:400).