

Anti-JAK1 (7B1) Mouse antibody



Catalog Number: 168126

Product name

WB, IP

Anti-JAK1 (7B1) Mouse antibody

Dilutions

Specificity

WB: 1/500

Human, Rat

Antibody description

Validations

JAK1 (7B1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Jak1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Clonality

Monoclonal


Ig Type

IgG1

Applications



Immunoprecipitation analysis of HeLa cell lysates using Jak1 mouse mAb (168126). Jak1 mouse mAb (200622-8B8) was used for the western blot



analysis (1:500 diluted)

Western blot analysis of extracts from PC12, C6, Ramos and Jurkat cell lysates using Jak1 mouse mAb (1:500 diluted). Predicted band size: 130kDa. Observed band size: 130kDa.