

Anti-Cystatin C (6A4) Mouse antibody



Catalog Number: 168159

Product name

Anti-Cystatin C (6A4) Mouse antibody

Specificity

Human

Antibody description

Cystatin C (6A4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant protein expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG1

Applications


WB, IHC-P, ELISA

Dilutions


WB: 1/1,000-2,000

IHC: 1/100-200 ELISA 1/10,000-20,000

Validations



Western blot detection of Cystatin C in NCI-H460, SKOV-3, PC3 cell lysates using Cystatin C (6A4) Mouse mAb(1:1000 diluted). Predicted band size: 14KDa. Observed band size: 14KDa.



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.