

Anti-RBP4 (8C12) Mouse antibody



Catalog Number: 168174

Product name

Monoclonal

Anti-RBP4 (8C12) Mouse antibody

Ig Type

Specificity

Mouse IgG1

Human

Applications

Antibody description

ELISA

RBP4 (8C12) Mouse monoclonal antibody

Validations

Preparation

Antigen: Purified recombinant human RBP4 protein fragments expressed in E.coli.


Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality



Standard Curve for RBP4: Capture Antibody Mouse mAb (8C12-C11) to RBP4 and Detector Antibody Mouse mAb(7C5-F6) to RBP4.