

# Anti-Lyn (1C8) Mouse antibody



Catalog Number: 168181

## Product name

Anti-Lyn (1C8) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human

## Applications

WB

## Antibody description

Lyn (1C8) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human Lyn protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Lyn in Ramos, Raji, K562 and HeLa cell lysates using Lyn mouse mAb(dilution 1:1000). Predicted band size: 56/53kDa. Observed band size: 56/53kDa.