

# Anti-beta Catenin (7E2) Mouse antibody



Catalog Number: 168202

## Product name

Anti-beta Catenin (7E2) Mouse antibody

WB: 1/1000-2000

IHC: 1/200-500

## Specificity

Human, Mouse, Rat

IF: 1/200

## Antibody description

beta Catenin (7E2) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of Catenin- $\beta$

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


IgG1

## Applications


WB, ICC/IF, IHC-P

## Dilutions

## Validations



Western blot analysis of 1) Hela, 2) 293T, 3) Mouse Liver Tissue, 4) Rat Liver Tissue, with  $\beta$ -Catenin Mouse mAb diluted at 1:2,000.



Western blot detection of  $\beta$ -Catenin in T47D, 3T3, COS7 and C6 cell lysates using  $\beta$ -Catenin mouse mAb (1:3000 diluted). Predicted band size: 92KDa. Observed band size: 92KDa.