

Anti-Transferrin (2D5) Mouse antibody



Catalog Number: 168217

Product name

Anti-Transferrin (2D5) Mouse antibody

Specificity

Human

Antibody description

Transferrin (2D5) Mouse monoclonal antibody

Preparation

Antigen: Transferrin human.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ELISA

Validations



Standard Curve for Transferrin: Capture Antibody Mouse mAb (2D5-E9) to Transferrin at 4u03bcg/ml and Detector Antibody Mouse mAb (4H5-H7) was used for detecting.



Observed Transferrin levels in several Human serum samples at different dilution.