

# Anti-Cystatin C (6F12) Mouse antibody



Catalog Number: 168225

## Product name

IgG1

Anti-Cystatin C (6F12) Mouse antibody

## Applications

## Specificity

sELISA Capture

Human

## Validations

## Antibody description

Cystatin C (6F12) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant human cystatin c protein.


## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Standard Curve for Human Cystatin C: Capture Antibody Mouse mAb [6F12-C7-D8] to Human Cystatin C at 4u03bcg/ml and Detector Antibody Mouse mAb [7F6-A5-F3]to Human Cystatin C at 0.5u03bcg/ml.