

Anti-Hsp27 (7E5) Mouse antibody



Catalog Number: 168231

Product name

Anti-Hsp27 (7E5) Mouse antibody

Specificity

Human, Monkey

Antibody description

Hsp27 (7E5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Hsp27 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/300

Validations



Western blot detection of Hsp27 in HeLa, HsHeLa, A431, COS7, Lncap and MCF7 cell lysates and using Hsp27 mouse mAb (1:1000 diluted). Predicted band size: 23KDa. Observed



band size:27KDa. Immunocytochemistry stain of HeLa using Hsp27 mouse mAb (1:300).