

Catalog Number: 168234

## Product name

Anti-PPP1A/PPP1CA (1C11) Mouse antibody

## Specificity

Human

## Antibody description

PPP1A/PPP1CA (1C11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PPP1A protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, IHC-P, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/300

## Validations




Immunocytochemistry stain of HeLa using PPP1A



mouse mAb (1:300).

Western blot detection of PPP1A in A431, Jurkat and HeLa whole cell lysates using PPP1A mouse mAb (1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.



IHC of paraffin-embedded human breast cancer using anti-PPP1A mouse mAb diluted 1/500-1/1000.