

# Anti-PPP1A/PPP1CA (1H1) Mouse antibody



Catalog Number: 168256

## Product name

Anti-PPP1A/PPP1CA (1H1) Mouse antibody

## Specificity

Human

## Antibody description

PPP1A/PPP1CA (1H1) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PPP1A protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/50

## Validations



Western blot detection of PPP1A in HepG2, MCF7 and HeLa cell lysates using PPP1A mouse mAb (1:1000 diluted). Predicted band size: 37KDa. Observed band size: 37KDa.



Immunocytochemistry of HeLa cells using anti-PPP1A mouse mAb diluted 1:50.