

# Anti-FABP-1 (6C6) Mouse antibody



Catalog Number: 168279

## Product name

Anti-FABP-1 (6C6) Mouse antibody

WB: 1/500 - 1/2000

ICC: 1/200 - 1/1000

## Specificity

Human, Mouse, Rat, Monkey

ELISA: 1/10000

## Antibody description

FABP-1 (6C6) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human GOT2 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions

## Validations



Western blot detection of GOT2 in MCF7, HCT116, Jurkat and Raji cell lysates using GOT2 mouse mAb (1:1000 diluted). Predicted band size: 47KDa. Observed band size: 47KDa.



Immunofluorescence analysis of PC-3 (left) and SK-BR-3 (right) cells using anti-GOT2 mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.