

Anti-NFI-C/CTF (4E10) Mouse antibody



Catalog Number: 168285

General Information

Product name:

Anti-NFI-C/CTF (4E10) Mouse antibody

Immunogen:

Applications:

WB, ICC/IF

Specificity:

Human

Antibody description:

NFI-C/CTF (4E10) Mouse monoclonal antibody

Properties

Preparation:

Antigen: Purified recombinant human NFIC protein fragments expressed in E.coli.

Formulation:

Ascites containing with 0.03% Proclin300 and 50% glycerol.

Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:

Monoclonal


Ig Type: IgG2b

Applications


WB: 1/200

;ICC: 1/200

Validations



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-NFIC mouse mAb (dilution 1:200).



Western blot detection of NFIC in HeLa NE cell lysates using NFIC mouse mAb (1:200 diluted). Predicted band size: 56 kDa. Observed band size: 56 kDa.