

# Anti-Ferritin Light Chain (1A2) Mouse antibody



Catalog Number: 168324

## Product name

Anti-Ferritin Light Chain (1A2) Mouse antibody

## Specificity

Human

## Antibody description

Ferritin Light Chain (1A2) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant full length of human ferritin light chain protein expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

IgG2b

## Applications

ELISA

## Validations



Indirect ELISA assay for anti-Ferritin Light Chain mouse mAb. Antigen coating concentration: 4ug/ml.