

Anti-DR5 (6F6) Mouse antibody



Catalog Number: 168335

Product name

Anti-DR5 (6F6) Mouse antibody

Dilutions

WB: 1/1000;Other 1/2000

Specificity

Human

Antibody description

DR5 (6F6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human DR5 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

IgG2b


Applications

WB

Validations



Western blot detection of DR5 in SW480, K562, HCT116 and HeLa cell lysates using DR5 mouse mAb (1:1000 diluted). Predicted band size: 40/48KDa. Observed band size: 40/48KDa.



Western blot detection of DR5 in MCF-7 and T47D cell lysates using DR5 mouse mAb (1:2000 diluted). Predicted band size: 40/48KDa. Observed band size: 40/48KDa.