

# Anti-GRP78 BiP (4F11) Mouse antibody



Catalog Number: 168342

## Product name

Anti-GRP78 BiP (4F11) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GRP78 BiP (4F11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human BiP/GRP78 protein expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG2b


## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000


## Validations



Immunoprecipitation of BiP/GRP78 from HeLa cell extracts using BiP/GRP78 Mouse mAb. Western blot was performed using BiP/GRP78 Mouse mAb.



Western blot analysis of extracts from 3T3, A549, C6 and HeLa cell lysates using BiP/GRP78 mouse mAb (1:1000 diluted). Predicted band size: 72KDa. Observed band size: 72KDa.



Immunofluorescent analysis of HeLa cells fixed with anhydrous methanol at -20°C and stained with BiP/GRP78 mouse mAb (dilution 1:200).