

Anti-Lck (2E5) Mouse antibody



Catalog Number: 168388

Product name

Anti-Lck (2E5) Mouse antibody

Specificity

Human

Antibody description

Lck (2E5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human Lck expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG2a

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot detection of Lck in Jurkat and Ramos cell lysates and using Lck mouse mAb (1:1000 diluted). Predicted band size: 56KDa. Observed band size: 56KDa.