

# Anti-p38 (6E9) Mouse antibody



Catalog Number: 168431

## Product name

WB: 1/500 - 1/2000

Anti-p38 (6E9) Mouse antibody

IHC: 1/200 - 1/1000

## Specificity

ELISA: 1/10000

Human, Mouse, Monkey, Rat

## Antibody description

p38 (6E9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human MAPK14 (AA: 299-360) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG1

## Applications


WB, IHC-P

## Dilutions

## Validations



Immunohistochemical analysis of paraffin-embedded endometrial cancer tissues using MAPK14 mouse mAb with DAB staining.



Western blot analysis of extracts from Jurkat, A549, MCF7, C6 and HeLa cell lysates using p38 mouse mAb (1:1000 diluted). Predicted band size: 41 kDa. Observed band size: 41 kDa.