

# Anti-beta 2 Microglobulin (4C3) Mouse antibody



Catalog Number: 168442

## Product name

Anti-beta 2 Microglobulin (4C3) Mouse antibody

## Specificity

Human

## Antibody description

beta 2 Microglobulin (4C3) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human  $\beta$ 2-microglobulin protein expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

IgG2b

## Applications

ELISA

## Validations



sELISA standard Curve for u03b22-microglobulin.