

# Anti-LIN28A (2C1) Mouse antibody



Catalog Number: 168454

## Product name

Anti-LIN28A (2C1) Mouse antibody

## Specificity

Human, Mouse

## Antibody description

LIN28A (2C1) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human LIN28A protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal in PBS(pH 7.4)containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG2a


## Applications

WB


## Dilutions

WB: 1/1000

## Validations



Western blot detection of LIN28 in NTERA and F9 cell lysates using LIN28A mouse mAb (1:1000 diluted).Predicted band size:26KDa.Observed



band size:26KDa.

Western blot detection of LIN28A in LIN28B and LIN28A recombinant antigen fragments and using LIN28A mouse mAb (1:1000 diluted).