

# Anti-Cardiac Troponin I (2D10) Mouse antibody



Catalog Number: 168455

## Product name

IgG2a

Anti-Cardiac Troponin I (2D10) Mouse antibody

## Applications

## Specificity

WB, ELISA

Human, Mouse, Rat

## Dilutions

## Antibody description

ELISA: 1/10000

Cardiac Troponin I (2D10) Mouse monoclonal antibody

;WB: 1/1000

## Preparation

## Validations

Antigen: Synthetic peptides corresponding to the sequence of human Cardiac Troponin I.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.


## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot analysis of extracts from Rat and mouse heart tissue using Troponin I (168455, dilution 1:1000) Mouse mAb. Predicted band size: 28kDa. Observed band size: 28kDa.