

Anti-ICAM1 (4F5) Mouse antibody



Catalog Number: 168471

Product name

Anti-ICAM1 (4F5) Mouse antibody

Dilutions

WB: 1/1000

Specificity

;ICC: 1/200

Human

Antibody description

ICAM1 (4F5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

IgG2b


Applications

WB, ICC/IF

Validations



Immunocytochemistry staining of SW480 cells fixed by anhydrous methanol for 2 h at -20°C and using anti-CD54(ICAM-1) mouse mAb



(dilution 1:200).

Western blot detection of CD54(ICAM-1) in SW480 and Raji cell lysates using CD54(ICAM-1) mouse mAb (1:1000 diluted). Predicted band size: 96KDa. Observed band size: 96KDa.