

# Anti-RICTOR Mouse antibody



Catalog Number: 168514

## Product name

Anti-RICTOR Mouse antibody

## Specificity

Human, Monkey, Mouse

## Antibody description

RICTOR Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human RICTOR expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB

## Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

## Validations



Western blot analysis using RICTOR mouse mAb against HeLa (1), PANC-1 (2), MOLT4 (3), and HepG2 (4) cell lysate.