

Anti-Hsc70 Mouse antibody



Catalog Number: 168515

Product name

Anti-Hsc70 Mouse antibody

Specificity

Human

Antibody description

Hsc70 Mouse monoclonal antibody

Preparation

Antigen: Recombinant Protein of Hsc 70.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 50% glycerol.

Storage

Store at -20°C

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of Hsc70 in HeLa,293 cell lysates using Hsc70 Mouse mAb(1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.