

Anti-Aquaporin 1 Rabbit antibody



Catalog Number: 168577

Product name

Anti-Aquaporin 1 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Aquaporin 1 Rabbit monoclonal antibody

Dilutions

1:1000-1:2000

Preparation


Antigen: A synthetic peptide of human Aquaporin 1

1:50-1:200

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Western blot detection of Aquaporin 1 in C6 cell lysates using Aquaporin 1 Rabbit mAb(1:1000 diluted).Predicted band size:29kDa.Observed band size:29kDa.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal