

Anti-Integrin alpha V Rabbit antibody



Catalog Number: 168586

Product name

Anti-Integrin alpha V Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Integrin alpha V Rabbit monoclonal antibody

Dilutions

1:1000

Preparation


1:20-1:500

Antigen: Recombinant protein of human Integrin alpha V

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Integrin alpha V in C6,3T3,HeLa cell lysates using Integrin alpha V Rabbit mAb(1:1000 diluted).Predicted band size:87kDa.Observed band size:135kDa.

Clonality

Monoclonal