

Anti-DNA Ligase IV Rabbit antibody



Catalog Number: 168607

Product name

Rabbit IgG

Anti-DNA Ligase IV Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF

Human

Dilutions

Antibody description

WB: 1/1000

DNA Ligase IV Rabbit monoclonal antibody

IHC: 1/20

Preparation

ICC/IF: 1/20

Antigen: Recombinant protein of human DNA Ligase IV

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of DNA Ligase IV in K562, U251MG, HepG2, MCF7 cell lysates using DNA Ligase IV Rabbit mAb(1:500 diluted). Predicted band size: 104kDa. Observed band size: 104kDa.