

Anti-ABAT/GABA-T Rabbit antibody



Catalog Number: 168644

Product name

Rabbit IgG

Anti-ABAT/GABA-T Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000

ABAT/GABA-T Rabbit monoclonal antibody

FC: 1/20

Preparation

IP: 1/20

Antigen: A synthetic peptide of human ABAT

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of ABAT/GABA-T in Rat Brain lysates using ABAT/GABA-T Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:56kDa.