

Anti-NAT1 Rabbit antibody



Catalog Number: 168683

Product name

Anti-NAT1 Rabbit antibody

Specificity

Human

Antibody description

NAT1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NAT1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000-1/5000

Validations



Western blot detection of NAT1 in C6,3T3,HeLa cell lysates using NAT1 Rabbit mAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.