

# Anti-p60 katanin Rabbit antibody



Catalog Number: 168687

## Product name

Anti-p60 katanin Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, ICC/IF

## Antibody description

p60 katanin Rabbit monoclonal antibody

## Dilutions

1:1000-1:5000

## Preparation


Antigen: A synthetic peptide of human p60 katanin

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

Western blot detection of p60 katanin in K562, HeLa cell lysates using p60 katanin Rabbit mAb(1:1000 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.