

# Anti-p53 Rabbit antibody



Catalog Number: 168699

## Product name

Anti-p53 Rabbit antibody

## Specificity

Human

## Antibody description

p53 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human p53

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000

ICC/IF: 1/50

## Validations



Western blot detection of p53 in Jurkat cell lysates using p53 Rabbit mAb(1:500 diluted).Predicted band size:44kDa.Observed band size:53kDa.