

Anti-DFNA5/GSDME Rabbit antibody



Catalog Number: 168715

Product name

Rabbit IgG

Anti-DFNA5/GSDME Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

DFNA5/GSDME Rabbit monoclonal antibody

FC: 1/20

Preparation


IP: 1/20

Antigen: Recombinant protein of human DFNA5/GSDME

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of DFNA5/GSDME in U251MG,U2OS,HepG2,Hela cell lysates using DFNA5/GSDME Rabbit mAb(1:500 diluted).Predicted band size:55kDa.Observed band size:55kDa.

Clonality

Monoclonal

Ig Type