

Anti-Ionotropic Glutamate receptor 2 Rabbit antibody



Catalog Number: 168734

Product name

Anti-Ionotropic Glutamate receptor 2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ionotropic Glutamate receptor 2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Ionotropic Glutamate receptor 2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

1:2000

1:50

1:20

Validations



Western blot detection of Ionotropic Glutamate receptor 2 in C6 cell lysates using Ionotropic Glutamate receptor 2 Rabbit mAb(1:1000 diluted).Predicted band size:99kDa.Observed band size:99kDa.