

Anti-Periostin Rabbit antibody



Catalog Number: 168742

Product name

Rabbit IgG

Anti-Periostin Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human

Dilutions

Antibody description

1:1000

Periostin Rabbit monoclonal antibody

1:200

Preparation


1:20

Antigen: A synthetic peptide of human Periostin

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of Periostin in C6 cell lysates using Periostin Rabbit mAb(1:1000 diluted).Predicted band size:93kDa.Observed band size:93kDa.

Ig Type