

# Anti-PAX8 Rabbit antibody



Catalog Number: 168743

## Product name

Rabbit IgG

Anti-PAX8 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

1:50

PAX8 Rabbit monoclonal antibody

1:50

## Preparation

1:20

Antigen: Recombinant protein of human PAX8

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of PAX8 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using PAX8 Rabbit mAb(1:1000 diluted). Predicted band size: 48 kDa; Observed MW: N/AkDa. Observed band size: N/AkDa.