

# Anti-Cystatin C Rabbit antibody



Catalog Number: 168755

## Product name

Anti-Cystatin C Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Cystatin C Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Cystatin C

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

1:2000-1:10000

1:50-1:200

## Validations



Western blot detection of Cystatin C in K562,3T3 cell lysates using Cystatin C Rabbit mAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:13kDa.