

Anti-COX2/cyclooxygenase2 Rabbit antibody



Catalog Number: 168786

Product name

Anti-COX2/cyclooxygenase2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, IHC-P

Antibody description

COX2/cyclooxygenase2 Rabbit monoclonal antibody

Dilutions

WB: 1/1000

IHC: 1/20

Preparation

Antigen: A synthetic peptide of human COX2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of COX-2/cyclooxygenase-2 in Rat Brain lysates using COX-2/cyclooxygenase-2 Rabbit mAb(1:500 diluted). Predicted band size:69kDa.Observed band size:69kDa.