

Anti-Vitamin D Binding protein Rabbit antibody



Catalog Number: 168813

Product name

Rabbit IgG

Anti-Vitamin D Binding protein Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC

Human

Dilutions

Antibody description

1:1000-1:5000

Vitamin D Binding protein Rabbit monoclonal antibody

1:20

Preparation

1:20

Antigen: A synthetic peptide of human Vitamin D Binding protein

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Vitamin D Binding protein in Rat Brain, 3T3, HeLa cell lysates using Vitamin D Binding protein Rabbit mAb (1:1000 diluted). Predicted band size: 53 kDa. Observed band size: 53 kDa.