

Anti-HMGB1 Rabbit antibody



Catalog Number: 168848

Product name

Rabbit IgG

Anti-HMGB1 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:2000

HMGB1 Rabbit monoclonal antibody

1:50

Preparation

1:50

Antigen: A synthetic peptide of human HMGB1

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of HMGB1 in K562,Rat Brain,C6,3T3,Hela cell lysates using HMGB1 Rabbit mAb(1:1000 diluted).Predicted band size:25kDa.Observed band size:25kDa.