

Anti-ATP6V1A Rabbit antibody



Catalog Number: 168918

Product name

Anti-ATP6V1A Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

1:1000

Antibody description

1:200

ATP6V1A Rabbit monoclonal antibody

1:50

Preparation

1:20

Antigen: Recombinant protein of human ATP6V1A

1:20

Formulation

Validations

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG



Western blot detection of ATP6V1A in K562, Rat Brain, 3T3 cell lysates using ATP6V1A Rabbit mAb(1:1000 diluted). Predicted band size: 68kDa. Observed band size: 68kDa.