

Anti-MCL1 Rabbit antibody



Catalog Number: 168962

Product name

Rabbit IgG

Anti-MCL1 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

MCL1 Rabbit monoclonal antibody

ICC/IF: 1/50

Preparation

FC: 1/50

Antigen: A synthetic peptide of human MCL1

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of MCL1 in C6 cell lysates using MCL1 Rabbit mAb(1:1000 diluted). Predicted band size: 37kDa. Observed band size: 40kDa.