

Anti-NR1D1 Rabbit antibody

Catalog Number: 168969

Product name

Anti-NR1D1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NR1D1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NR1D1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000


ICC/IF: 1/20

FC: 1/20

Validations



Immunocytochemistry of NR1D1(green) in HeLa cells using NR1D1 Rabbit mAb at dilution 1/50, and DAPI(blue)



and DAPI(blue)

Western blot detection of NR1D1 in K562,C6,3T3,HeLa cell lysates using NR1D1 Rabbit mAb(1:1000 diluted).Predicted band size:67kDa.Observed band size:67kDa.