

# Anti-ZNF364/RFN115 Rabbit antibody



Catalog Number: 169028

## Product name

Anti-ZNF364/RFN115 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

ZNF364/RFN115 Rabbit monoclonal antibody

## Dilutions

1:1000


## Preparation

Antigen: A synthetic peptide of human ZNF364

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of ZNF364/RFN115 in Jurkat, HeLa cell lysates using ZNF364/RFN115 Rabbit mAb(1:1000 diluted). Predicted band size:34kDa. Observed band size:34kDa.

## Clonality

Monoclonal